

Bullet Protocol

Preparation of Gene Gun Bullets and Biolistic Transfection of Neurons in Slice Culture
<http://www.jove.com/video/675/preparation-of-gene-gun-bullets-and-biolistic-transfection-of-neurons-in-slice-culture>

Au solution (200mg/mL in 50% glycerol, size 1.6 μ m Gold Microcarriers - 165-2264)

Spermidine (50mM stock in deep freezer)

1M CaCl₂ (help to stick DNA to Au)

EtOH, anhydrous (Sigma 459836). Use 10mL syringe.

20mg/mL PVP in dw (make fresh, sonication will help to dissolve). PVP will help the Au to stick to the tube. *polyvinylpyrrolidone*

1. Prepare DNA solution. (about 3 μ g / construct). Total DNA volume should be less than 6 μ L. DNA should be maxi prep. If using two different construct, mix well
2. Combine 15 μ L of Au solution with 17 μ L of 50mM spermidine. Au solution should be well vortexed and sonicated. After combined, Sonicate 5sec, 3times. *20 (24) 22 (26) cut end of tip -*
3. Combine DNA and Au/spermidine solution.
4. Drop wise 17 μ L of 1M CaCl₂ while vortexing.
5. Incubate for 3min, RT *20*
6. Cfg 6,000rpm, 30sec and remove sup *absol. - use syringe*
7. Wash and spin (full speed for 30sec) three times with 1mL of 100% EtOH to remove all traces of water. Sonicate for 1sec (don't over Sonicate) before each centrifugation making sure there is no aggregated Au prior to spinning. *- brief*
8. In a separate tube prepare 0.05mg/mL PVP EtOH by adding 3.75 μ L of 20mg/mL PVP (dessicated) in 1.5mL EtOH. *7.5*
9. Combine 1.5mL 0.05mg/mL PVP with Au pellet. (use less than 1.5mL may concentrate the Au particle per tube) *→ make better stick on tubing*
10. Vortex immediately prior to loading into tube.
11. Use syringe to slowly suck Au into two foot piece of tubing (1.5mL will fill about half of the tube; mark the Au boundary on the tube. Leave a few inches dry from the end). Careful to avoid bubbles.
12. Let Au settle to bottom of tube for 10min